

Catalogue of Autoclave Bags



Autoclavable Bags

Products features:

Superior strength HDEP bag, are free from chlorine.

Size: 457x457, 508x610, 610x660, 660x762, 762x914, and 1016x1270 in millimeter (customize as per requirement)

Thickness: 55 to 75 micron (as per guideline)

Color: Red, White, Transparent.

Not deform in autoclave machine.

May be imprinted with the biohazard symbol

Environment Friendly

Can be steam sterilization at 121 °C- 15 minutes

Decontamination and Disposal of contaminated material made easy.

Transparent – easy to see through.

Available in multiple packing of 50 bags.

These bags are Autoclavable bag used for waste disposal and recommended for disposal of pathological/ clinical or contaminated material. It can also be used for sterilization of glass ware or plastic wares.

Application:

Used in hospitals, clinical laboratories, microbiological laboratories, Biotechnology laboratories, molecular biology laboratories etc. for disposal of clinical and pathological or contaminated material.

Direction for use:

1. Place all the clinical / contaminated material in the bag.
2. Seal the bag.
3. Autoclave at 121°C for 15 minutes.
4. Discard the autoclaved bag with decontaminated material as per local law